

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat TLE4 Polyclonal Antibody

多克隆抗体

产品货号: PA06034

---

第三版

**描述:**TLE4 (Transducin-like enhancer protein 4) is a member of the Groucho family of transcriptional corepressors and inhibits the transcriptional activation mediated by Pax5 and by beta-catenin and TCF family members in Wnt signaling.

**抗原:**Synthesized peptide derived from the Internal region of human TLE4

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**In all tissues examined, mostly in brain, and muscle.