

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TOLLIP Polyclonal Antibody

多克隆抗体

产品货号: PA06043

第三版

**描述:** This gene encodes a ubiquitin-binding protein that interacts with several Toll-like receptor (TLR) signaling cascade components. The encoded protein regulates inflammatory signaling and is involved in interleukin-1 receptor trafficking and in the turnover of IL1R-associated kinase. Several transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthesized peptide derived from the Internal region of human Tollip.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**