

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TOMM20 Polyclonal Antibody

多克隆抗体

产品货号: PA06044

第三版

描述:TOMM20 (Translocase Of Outer Mitochondrial Membrane 20) is a Protein Coding gene. Diseases associated with TOMM20 include Perry Syndrome and Infantile Neuroaxonal Dystrophy 1. Among its related pathways are Deubiquitination and Neuroscience. GO annotations related to this gene include unfolded protein binding and P-P-bond-hydrolysis-driven protein transmembrane transporter activity. An important paralog of this gene is TOMM20L.

抗原:Synthesized peptide derived from the Internal region of human Tom20

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: