

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human UBA2 Polyclonal Antibody

多克隆抗体

产品货号: PA06081

第三版

描述:The heterodimer acts as a E1 ligase for SUMO1, SUMO2, SUMO3, and probably SUMO4. It mediates ATP-dependent activation of SUMO proteins followed by formation of a thioester bond between a SUMO protein and a conserved active site cysteine residue on UBA2/SAE2.

抗原:Synthesized peptide derived from the C-terminal region of human UBA2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: