

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human UBE3C Polyclonal Antibody

多克隆抗体

产品货号: PA06087

第三版

描述:UBE3C (Ubiquitin Protein Ligase E3C) is a Protein Coding gene. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. GO annotations related to this gene include ligase activity and ubiquitin-protein transferase activity. An important paralog of this gene is UBE3B.

抗原:Synthesized peptide derived from the Internal region of human UBE3C.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in skeletal muscle. Detected at much lower levels in kidney and pancreas.