

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human USP38 Polyclonal Antibody

多克隆抗体

产品货号: PA06101

第三版

**描述:** USP38 (ubiquitin specific peptidase 38), also known as HP43.8KD, is a 1,042 amino acid deubiquitinating enzyme that participates in the Ub pathway, exhibiting a preference towards 'Lys-63'-linked Ubiquitin chains. A member of the peptidase C19 family, the catalytic activity of USP38 involves a combination of the ubiquitin carboxyl-terminal thiolester and water to produce ubiquitin and a thiol. USP38 is highly expressed in skeletal muscle and adrenal gland.

**抗原:** Synthesized peptide derived from the Internal region of human USP38

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:**