

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human USP44 Polyclonal Antibody

多克隆抗体

产品货号: PA06105

第三版

描述:Deubiquitinase that plays a key role in the spindle checkpoint by preventing premature anaphase onset. Acts by specifically mediating deubiquitination of CDC20, a negative regulator of the anaphase promoting complex/cyclosome (APC/C). Deubiquitination of CDC20 leads to stabilize the MAD2L1-CDC20-APC/C ternary complex (also named mitotic checkpoint), thereby preventing premature activation of the APC/C.

抗原:Synthesized peptide derived from the Internal region of human USP44

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in testis.