

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat VAV1 Polyclonal Antibody

多克隆抗体

产品货号: PA06116

第三版

描述: This gene is a member of the VAV gene family. The VAV proteins are guanine nucleotide exchange factors (GEFs) for Rho family GTPases that activate pathways leading to actin cytoskeletal rearrangements and transcriptional alterations. The encoded protein is important in hematopoiesis, playing a role in T-cell and B-cell development and activation. The encoded protein has been identified as the specific binding partner of Nef proteins from HIV-1. Coexpression and binding of these partners initiates profound morphological changes, cytoskeletal rearrangements and the JNK/SAPK signaling cascade, leading to increased levels of viral transcription and replication. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. VAV1 (Vav Guanine Nucleotide Exchange Factor 1) is a Protein Coding gene. Diseases associated with VAV1 include Spherocytosis, Type 1 and Common Variable Immunodeficiency. Among its related pathways are ICos-ICosL Pathway in T-Helper Cell and Type I Interferon Signaling Pathways. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and Rho guanyl-nucleotide exchange factor activity. An important paralog of this gene is VAV3.

抗原: Synthesized peptide derived from human Vav1 around the non-phosphorylation site of Tyr174.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed in hematopoietic cells but not in other cell types.