

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human WDHD1 Polyclonal Antibody

多克隆抗体

产品货号: PA06136

第三版

描述:The protein encoded by this gene contains multiple N-terminal WD40 domains and a C-terminal high mobility group (HMG) box. WD40 domains are found in a variety of eukaryotic proteins and may function as adaptor/regulatory modules in signal transduction, pre-mRNA processing and cytoskeleton assembly. HMG boxes are found in many eukaryotic proteins involved in chromatin assembly, transcription and replication. Alternative splicing results in two transcript variants encoding different isoforms.WDHD1 (WD Repeat And HMG-Box DNA Binding Protein 1) is a Protein Coding gene. Diseases associated with WDHD1 include Trench Fever. Among its related pathways are Mesodermal Commitment Pathway. GO annotations related to this gene include RNA binding and RNA polymerase binding.

抗原:Synthesized peptide derived from the Internal region of human WDHD1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: