

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human XRCC3 Polyclonal Antibody

多克隆抗体

产品货号: PA06147

第三版

描述:Involved in the homologous recombination repair (HRR) pathway of double-stranded DNA, thought to repair chromosomal fragmentation, translocations and deletions. Plays a role in regulating mitochondrial DNA copy number under conditions of oxidative stress in the presence of RAD51 and RAD51C.

抗原:Synthesized peptide derived from the Internal region of human XRCC3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: