

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human YEATS2 Polyclonal Antibody

多克隆抗体

产品货号: PA06152

第三版

描述:The protein encoded by this gene is a scaffolding subunit of the ATAC complex, which is a complex with acetyltransferase activity on histones H3 and H4. Alternative splicing results in multiple transcript variants encoding different isoforms.YEATS2 (YEATS Domain Containing 2) is a Protein Coding gene. Among its related pathways are Chromatin organization. GO annotations related to this gene include TBP-class protein binding. An important paralog of this gene is YEATS4.

抗原:Synthesized peptide derived from the Internal region of human YEATS2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: