

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZAR1 Polyclonal Antibody

多克隆抗体

产品货号: PA06156

第三版

描述:This maternal effect gene is oocyte-specific and encodes a protein that is thought to function in the initiation of embryogenesis. A similar protein in mouse is required for female fertility. ZAR1 (Zygote Arrest 1) is a Protein Coding gene. Among its related pathways are Preimplantation Embryo. An important paralog of this gene is ZAR1L.

抗原:Synthesized peptide derived from the Internal region of human ZAR1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ovary and testis.