

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey ZC3H11A Polyclonal Antibody

多克隆抗体

产品货号: PA06159

第三版

描述:ZC3H11A (Zinc Finger CCCH-Type Containing 11A) is a Protein Coding gene. Among its related pathways are Gene Expression and Transport of Mature Transcript to Cytoplasm. GO annotations related to this gene include poly(A) RNA binding. An important paralog of this gene is ZC3H11B.

抗原:Synthesized peptide derived from the Internal region of human ZC3H11A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: