

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC39A4 Polyclonal Antibody

多克隆抗体

产品货号: PA06163

第三版

描述: This gene encodes a member of the zinc/iron-regulated transporter-like protein (ZIP) family. The encoded protein localizes to cell membranes and is required for zinc uptake in the intestine. Mutations in this gene result in acrodermatitis enteropathica. Multiple transcript variants encoding different isoforms have been found for this gene. SLC39A4 (Solute Carrier Family 39 Member 4) is a Protein Coding gene. Diseases associated with SLC39A4 include Acrodermatitis Enteropathica and Acrodermatitis. Among its related pathways are Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds and Senescence and Autophagy in Cancer. GO annotations related to this gene include metal ion transmembrane transporter activity and zinc ion transmembrane transporter activity. An important paralog of this gene is SLC39A12.

抗原: Synthesized peptide derived from the Internal region of human ZIP4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in kidney, small intestine, stomach, colon, jejunum and duodenum.