
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ZnT-9 Polyclonal Antibody

多克隆抗体

产品货号: PA06187

第三版

描述:SLC30A9 (Solute Carrier Family 30 Member 9) is a Protein Coding gene. Among its related pathways are Metal ion SLC transporters. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and ligand-dependent nuclear receptor transcription coactivator activity.

抗原:Synthesized peptide derived from the Internal region of human ZnT-9

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed in fetal and adult tissues and cancer cell lines.