

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZP1 Polyclonal Antibody

多克隆抗体

产品货号: PA06189

第三版

**描述:**The zona pellucida is an extracellular matrix that surrounds the oocyte and early embryo. It is composed primarily of three or four glycoproteins with various functions during fertilization and preimplantation development. The protein encoded by this gene ensures the structural integrity of the zona pellucida. Mutations in this gene are a cause of oocyte maturation defect and infertility.

**抗原:**Synthesized peptide derived from the Internal region of human ZP1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Oocytes.