

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CHAF1A Polyclonal Antibody

多克隆抗体

产品货号: PA06193

第三版

**描述:**Chromatin assembly factor I (CAF1) is a nuclear complex consisting of p50, p60 (CHAF1B; MIM 601245), and p150 (CHAF1A) subunits that assembles histone octamers onto replicating DNA in vitro. CHAF1A (Chromatin Assembly Factor 1 Subunit A) is a Protein Coding gene. Among its related pathways are Autophagy Pathway and Chromatin Regulation / Acetylation. GO annotations related to this gene include identical protein binding and unfolded protein binding.

**抗原:**Synthesized peptide derived from the Internal region of human CAF-1 p150

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**