

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat AAAS Polyclonal Antibody

多克隆抗体

产品货号: PA06211

第三版

描述: Adracalin (AAAS) is expressed in both neuroendocrine and cerebral structures and may function in the normal development of the peripheral and central nervous system. It localizes to nuclear pore complexes (NPCs), large multiprotein assemblies that are the sole sites of nucleocytoplasmic transport. Defects in AAAS are the cause of achalasia-addisonianism-alacrima syndrome (AAA syndrome); also known as triple-A syndrome or Allgrove syndrome.

抗原: Synthesized peptide derived from the C-terminal region of human Aladin

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed. Particularly abundant expression is found in cerebellum, corpus callosum, adrenal gland, pituitary gland, gastrointestinal structures and fetal lung.