

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human A1BG Polyclonal Antibody

多克隆抗体

产品货号: PA06222

第三版

描述:Alpha-1-B glycoprotein is a protein in humans that is encoded by the A1BG gene. The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins.

抗原:Synthesized peptide derived from the Internal region of human A1BG

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Plasma.