

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SV2C Polyclonal Antibody

多克隆抗体

产品货号: PA06239

第三版

**描述:**SV2C (Synaptic Vesicle Glycoprotein 2C) is a Protein Coding gene. Diseases associated with SV2C include Alcohol-Related Birth Defect and Fetal Alcohol Spectrum Disorder. Among its related pathways are Uptake and actions of bacterial toxins and HIV Life Cycle. GO annotations related to this gene include transporter activity and transmembrane transporter activity. An important paralog of this gene is SV2A.

**抗原:**Synthesized peptide derived from the Internal region of human SV2C.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**