

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MPHOSPH10 Polyclonal Antibody

多克隆抗体

产品货号: PA06241

第三版

描述: This gene encodes a protein that is phosphorylated during mitosis. The protein localizes to the nucleolus during interphase and to the chromosomes during M phase. The protein associates with the U3 small nucleolar ribonucleoprotein 60-80S complexes and may be involved in pre-rRNA processing. MPHOSPH10 (M-Phase Phosphoprotein 10) is a Protein Coding gene. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. GO annotations related to this gene include poly(A) RNA binding.

抗原: Synthesized peptide derived from the Internal region of human MPP10

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: