

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STRN Polyclonal Antibody

多克隆抗体

产品货号: PA06250

第三版

描述:STRN (Striatin) is a Protein Coding gene. Diseases associated with STRN include Cherubism and Cerebral Cavernous Malformations 3. Among its related pathways are Plasma membrane estrogen receptor signaling and Neurophysiological process Glutamate regulation of Dopamine D1A receptor signaling. GO annotations related to this gene include protein complex binding and estrogen receptor binding. An important paralog of this gene is STRN3.

抗原:Synthesized peptide derived from the C-terminal region of human Striatin.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Preferentially expressed in brain.