

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SLC6A17 Polyclonal Antibody

多克隆抗体

产品货号: PA06255

第三版

**描述:**The protein encoded by this gene is a member of the SLC6 family of transporters, which are responsible for the presynaptic uptake of most neurotransmitters. The encoded vesicular transporter is selective for proline, glycine, leucine and alanine. Defects in this gene are a cause of autosomal recessive mental retardation-48.SLC6A17 (Solute Carrier Family 6 Member 17) is a Protein Coding gene. Diseases associated with SLC6A17 include Mental Retardation, Autosomal Recessive 48. Among its related pathways are NRF2 pathway. GO annotations related to this gene include neurotransmitter:sodium symporter activity. An important paralog of this gene is SLC6A15.

**抗原:**Synthesized peptide derived from the Internal region of human NTT4

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**