

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PAOX Polyclonal Antibody

多克隆抗体

产品货号: PA06256

第三版

**描述:**PAOX (Polyamine Oxidase) is a Protein Coding gene. Among its related pathways are CDK-mediated phosphorylation and removal of Cdc6 and Peroxisome. GO annotations related to this gene include receptor binding and polyamine oxidase activity. An important paralog of this gene is SMOX.

**抗原:**Synthesized peptide derived from the Internal region of human PAOX

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Widely expressed. Not detected in spleen. Expressed at lower level in neoplastic tissues.