

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SCML2 Polyclonal Antibody

多克隆抗体

产品货号: PA06259

第三版

描述: This gene encodes a member of the Polycomb group proteins. These proteins form the Polycomb repressive complexes which are involved in transcriptional repression. The encoded protein binds histone peptides that are monomethylated at lysine residues and may be involved in regulating homeotic gene expression during development. SCML2 (Sex Comb On Midleg Like 2 (Drosophila)) is a Protein Coding gene. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and sequence-specific DNA binding. An important paralog of this gene is SCMH1.

抗原: Synthesized peptide derived from the Internal region of human SCML2.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in placenta, thymus and testis. Detected at lower levels in brain, liver, skeletal muscle, pancreas and ovary.