

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat WNT16 Polyclonal Antibody

多克隆抗体

产品货号: PA06263

第三版

**描述:**The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It contains two transcript variants diverging at the 5' termini. These two variants are proposed to be the products of separate promoters and not to be splice variants from a single promoter. They are differentially expressed in normal tissues, one of which (variant 2) is expressed at significant levels only in the pancreas, whereas another one (variant 1) is expressed more ubiquitously with highest levels in adult kidney, placenta, brain, heart, and spleen. WNT16 (Wnt Family Member 16) is a Protein Coding gene. Diseases associated with WNT16 include Nodular Basal Cell Carcinoma and Basal Cell Carcinoma. Among its related pathways are Wnt Signaling Pathways: beta-Catenin-independent Wnt/Ca<sup>2+</sup> Signaling and Other Non-canonical Wnt Signaling Pathways and Transcription Androgen Receptor nuclear signaling. GO annotations related to this gene include receptor binding and frizzled binding. An important paralog of this gene is WNT7B.

**抗原:**Synthesized peptide derived from the Internal region of human Wnt-16.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Isoform Wnt-16b is expressed in peripheral lymphoid organs such as spleen, appendix, and lymph nodes, in kidney but not in bone marrow. Isoform Wnt-16a is expressed at significant levels only in the pancreas.