

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat RRP7A Polyclonal Antibody

多克隆抗体

产品货号: PA06268

---

第三版

**描述:**RRP7A (Ribosomal RNA Processing 7 Homolog A) is a Protein Coding gene. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. GO annotations related to this gene include poly(A) RNA binding and nucleotide binding.

**抗原:**Synthesized peptide derived from the C-terminal region of human RRP7A.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**