

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HSP77/76 Polyclonal Antibody

多克隆抗体

产品货号: PA06269

第三版

描述:HSPA7 (Heat Shock Protein Family A (Hsp70) Member 7) is a Pseudogene. Among its related pathways are Mechanisms of CFTR activation by S-nitrosoglutathione (normal and CF) and Proteolysis_Putative ubiquitin pathway.

抗原:Synthesized peptide derived from the Internal region of human HSP77/76.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: