

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SYCP2 Polyclonal Antibody

多克隆抗体

产品货号: PA06271

第三版

描述:The synaptonemal complex is a proteinaceous structure that links homologous chromosomes during the prophase of meiosis. The protein encoded by this gene is a major component of the synaptonemal complex and may bind DNA at scaffold attachment regions. The encoded protein requires synaptonemal complex protein 3, but not 1, for inclusion in the synaptonemal complex.

抗原:Synthesized peptide derived from the Internal region of human SCP-2.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: