

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TXNDC17 Polyclonal Antibody

多克隆抗体

产品货号: PA06285

第三版

描述:Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF-alpha signaling and NF-kappa-B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.

抗原:Synthesized peptide derived from the Internal region of human TRP14

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed in cell lines.