

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PREP Polyclonal Antibody

多克隆抗体

产品货号: PA06286

第三版

描述:The protein encoded by this gene is a cytosolic prolyl endopeptidase that cleaves peptide bonds on the C-terminal side of prolyl residues within peptides that are up to approximately 30 amino acids long. Prolyl endopeptidases have been reported to be involved in the maturation and degradation of peptide hormones and neuropeptides. PREP (Prolyl Endopeptidase) is a Protein Coding gene. Diseases associated with PREP include Celiac Disease and Post-Traumatic Stress Disorder. Among its related pathways are Agents Acting on the Renin-Angiotensin System Pathway, Pharmacodynamics. GO annotations related to this gene include serine-type endopeptidase activity and serine-type exopeptidase activity. An important paralog of this gene is PREPL.

抗原:Synthesized peptide derived from the N-terminal region of human PREP

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: