

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL4 Polyclonal Antibody

多克隆抗体

产品货号: PA06298

第三版

描述:Participates in at least several B-cell activation processes as well as of other cell types. It is a costimulator of DNA-synthesis. It induces the expression of class II MHC molecules on resting B-cells. It enhances both secretion and cell surface expression of IgE and IgG1. It also regulates the expression of the low affinity Fc receptor for IgE (CD23) on both lymphocytes and monocytes.

抗原:Synthesized peptide derived from the Internal region of human IL-4.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: