

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EDN1 Polyclonal Antibody

多克隆抗体

产品货号: PA06300

第三版

描述:The protein encoded by this gene is proteolytically processed to release a secreted peptide termed endothelin 1. This peptide is a potent vasoconstrictor and is produced by vascular endothelial cells. Endothelin 1 also can affect the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.

抗原:Synthesized peptide derived from the Internal region of human ET-1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in lung, placental stem villi vessels and in cultured placental vascular smooth muscle cells.