

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat FGF12 Polyclonal Antibody

多克隆抗体

产品货号: PA06312

第三版

描述:Involved in nervous system development and function. Promote neuronal excitability by elevating the voltage dependence of neuronal sodium channel SCN8A fast inactivation

抗原:Synthesized peptide derived from the Internal region of human FGF-12.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Brain, eye and testis; highly expressed in embryonic retina, olfactory epithelium, olfactory bulb, and in a segmental pattern of the body wall; in adult olfactory bulb, less in cerebellum, deep cerebellar nuclei, cortex and multiple midbrain structures.