

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HSP40-4 Polyclonal Antibody

多克隆抗体

产品货号: PA06313

第三版

**描述:**Co-chaperone of Hsc70. Seems to play a role in protein import into mitochondria.

**抗原:**Synthesized peptide derived from the Internal region of human HSP40-4.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Ubiquitous. Isoform 2 is highly expressed in testis and lung, but detected at low levels in thymus, prostate, colon and liver.