

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CPE Polyclonal Antibody

多克隆抗体

产品货号: PA06317

第三版

描述:Removes residual C-terminal Arg or Lys remaining after initial endoprotease cleavage during prohormone processing. Processes proinsulin.

抗原:Synthesized peptide derived from the C-terminal region of human CPE.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Isoform 2, but not isoform 1, is overexpressed in hepatocellular carcinoma (at protein level), as well as in other tumors, including pheochromocytomas and paragangliomas.