

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LCN1 Polyclonal Antibody

多克隆抗体

产品货号: PA06366

第三版

描述: This gene encodes a member of the lipocalin family of small secretory proteins. Lipocalins are extracellular transport proteins that bind to a variety of hydrophobic ligands. The encoded protein is the primary lipid binding protein in tears and is overproduced in response to multiple stimuli including infection and stress. The encoded protein may be a marker for chromosome aneuploidy as well as an autoantigen in Sjogren's syndrome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and two pseudogenes of this gene are also located on the long arm of chromosome 9.

抗原: Synthesized peptide derived from the Internal region of human Lipocalin-1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Mainly expressed in lachrymal and salivary glands. Also expressed in the prostate.