

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PROC Polyclonal Antibody

多克隆抗体

产品货号: PA06368

第三版

**描述:** This gene encodes a vitamin K-dependent plasma glycoprotein. The encoded protein is cleaved to its activated form by the thrombin-thrombomodulin complex. This activated form contains a serine protease domain and functions in degradation of the activated forms of coagulation factors V and VIII. Mutations in this gene have been associated with thrombophilia due to protein C deficiency, neonatal purpura fulminans, and recurrent venous thrombosis. [provided by RefSeq, Dec 2009]

**抗原:** Synthesized peptide derived from the Internal region of human Protein C.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Plasma; synthesized in the liver.