

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD1E Polyclonal Antibody

多克隆抗体

产品货号: PA06387

第三版

**描述:**This gene encodes a member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes within Golgi compartments, endosomes, and lysosomes, and is cleaved into a stable soluble form. The soluble form is required for the intracellular processing of some glycolipids into a form that can be presented by other CD1 family members. Many alternatively spliced transcript variants encoding different isoforms have been described. Additional transcript variants have been found; however, their biological validity has not been determined. CD1E (CD1e Molecule) is a Protein Coding gene. Diseases associated with CD1E include Mycobacterium Malmoense and Autoimmune Disease Of Central Nervous System. Among its related pathways are Tight junction and Hematopoietic cell lineage. GO annotations related to this gene include lipid binding and lipopeptide binding. An important paralog of this gene is CD1B.

**抗原:**Synthesized peptide derived from the C-terminal region of human CD1E.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed on cortical thymocytes, dendritic cells, Langerhans cells, on certain T-cell leukemias, and in various other tissues.