

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD248 Polyclonal Antibody

多克隆抗体

产品货号: PA06409

第三版

描述:CD248 (CD248 Molecule) is a Protein Coding gene. GO annotations related to this gene include calcium ion binding and extracellular matrix binding. An important paralog of this gene is CD93.

抗原:Synthesized peptide derived from the Internal region of human CD248.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in tumor endothelial cells but absent or barely detectable in normal endothelial cells. Expressed in metastatic lesions of the liver and during angiogenesis of corpus luteum formation and wound healing. Expressed in vascular endothelial cells of malignant tumors but not in normal blood vessels. Expressed in stromal fibroblasts.