

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD229 Polyclonal Antibody

多克隆抗体

产品货号: PA06413

第三版

描述:LY9 belongs to the SLAM family of immunomodulatory receptors (see SLAMF1; MIM 603492) and interacts with the adaptor molecule SAP.LY9 (Lymphocyte Antigen 9) is a Protein Coding gene. Diseases associated with LY9 include Systemic Lupus Erythematosus. Among its related pathways are Hematopoietic Stem Cell Differentiation Pathways and Lineage-specific Markers. GO annotations related to this gene include receptor activity. An important paralog of this gene is CD84.

抗原:Synthesized peptide derived from the Internal region of human CD229.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Increased surface expression on T-cells of systemic lupus erythematosus (SLE) patients.