

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GDF5 Polyclonal Antibody

多克隆抗体

产品货号: PA06440

第三版

描述: This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This protein regulates the development of numerous tissue and cell types, including cartilage, joints, brown fat, teeth, and the growth of neuronal axons and dendrites. Mutations in this gene are associated with acromesomelic dysplasia, brachydactyly, chondrodysplasia, multiple synostoses syndrome, proximal symphalangism, and susceptibility to osteoarthritis. GDF5 (Growth Differentiation Factor 5) is a Protein Coding gene. Diseases associated with GDF5 include Du Pan Syndrome and Chondrodysplasia, Grebe Type. Among its related pathways are GPCR Pathway and Phospholipase-C Pathway. GO annotations related to this gene include identical protein binding and cytokine activity. An important paralog of this gene is GDF6.

抗原: Synthesized peptide derived from the Internal region of human GDF-5.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Predominantly expressed in long bones during embryonic development.