

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CMTM6 Polyclonal Antibody

多克隆抗体

产品货号: PA06445

第三版

**描述:** This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown.

**抗原:** Synthesized peptide derived from the C-terminal region of human CMTM6.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Expressed in the leukocytes, placenta and testis.