

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINB3 Polyclonal Antibody

多克隆抗体

产品货号: PA06447

第三版

描述:SERPINB3(serpin peptidase inhibitor, clade B (ovalbumin), member 3), also known as SCCA1, increase in advancing stages of some cervical, lung, esophageal and head and neck cancers). Human chromosome position 18q21.3 contains a cluster of serpins, including a tandem duplication of the SCCA gene, plasminogen activator inhibitor type 2, and maspin. SCCA is transcribed by two nearly identical genes (SCCA1 and SCCA2), and is mainly produced as SCCA1. The human SCCA1 gene encodes a 390 amino acid protein that was originally isolated from a metastatic cervical squamous cell carcinoma.

抗原:Synthesized peptide derived from the Internal region of human SCCA1/2.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Squamous cells. Expressed in some hepatocellular carcinoma (at protein level).