

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZKSCAN4 Polyclonal Antibody

多克隆抗体

产品货号: PA06452

第三版

**描述:**ZKSCAN4 (Zinc Finger With KRAB And SCAN Domains 4) is a Protein Coding gene. Among its related pathways are Gene Expression. GO annotations related to this gene include nucleic acid binding and transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is ZKSCAN3.

**抗原:**Synthesized peptide derived from the Internal region of human ZKSCAN4.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in adult heart, brain, placenta, lung and kidney, but not in adult liver and skeletal muscle. In 17-day old embryo, detected in liver, skeletal muscle, brain, heart and small intestine.