

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CK-8 Polyclonal Antibody

多克隆抗体

产品货号: PA06462

第三版

描述: Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

抗原: Synthesized peptide derived from human Cytokeratin 8 around the non-acetylation site of K483.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Observed in muscle fibers accumulating in the costameres of myoplasm at the sarcolemma membrane in structures that contain dystrophin and spectrin. Expressed in gingival mucosa and hard palate of the oral cavity.