

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CSF1 Polyclonal Antibody

多克隆抗体

产品货号: PA06489

第三版

描述:The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants.

抗原:Synthesized peptide derived from the C-terminal region of human M-CSF.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: