

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human STRA8 Polyclonal Antibody

多克隆抗体

产品货号: PA06492

第三版

**描述:** This gene encodes a retinoic acid-responsive protein. A homologous protein in mouse has been shown to be involved in the regulation of meiotic initiation in both spermatogenesis and oogenesis, though feature differences between the mouse and human proteins suggest that these homologs are not entirely functionally equivalent. It is thought that this gene may play a role in spermatogenesis in humans.

**抗原:** Synthesized peptide derived from the Internal region of human Stra8.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed specifically in testis.