

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINA1 Polyclonal Antibody

多克隆抗体

产品货号: PA06494

第三版

描述:The protein encoded by this gene is secreted and is a serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can cause emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene.

抗原:Synthesized peptide derived from the N-terminal region of human AAT.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous. Expressed in leukocytes and plasma.