

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PPP1R15B Polyclonal Antibody

多克隆抗体

产品货号: PA06502

第三版

描述: This gene encodes a protein phosphatase I-interacting protein that promotes the dephosphorylation of eukaryotic translation initiation factor 2A to regulate translation under conditions of cellular stress. The transcribed messenger RNA contains two upstream open reading frames (ORFs) that repress translation of the main protein coding ORF under normal conditions, while the protein coding ORF is expressed at high levels in response to stress. Continual translation of the mRNA under conditions of eukaryotic translation initiation factor 2A inactivation is thought to create a feedback loop for reactivation of the gene during recovery from stress. In addition, it has been shown that this protein plays a role in membrane traffic that is independent of translation and that it is required for exocytosis from erythroleukemia cells. Allelic variants of this gene are associated with microcephaly, short stature, and impaired glucose metabolism. PPP1R15B (Protein Phosphatase 1 Regulatory Subunit 15B) is a Protein Coding gene. Diseases associated with PPP1R15B include Microcephaly, Short Stature, And Impaired Glucose Metabolism 2 and Primary Microcephaly-Mild Intellectual Disability-Young-Onset Diabetes Syndrome. GO annotations related to this gene include protein serine/threonine phosphatase activity.

抗原: Synthesized peptide derived from the Internal region of human PPP1R15B.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: